

Pen-type electronic pipette with smartphone APP
Pen-type electronic pipette with wireless communication function
pipetty Smart APP (Dedicated application)



Operation Manual

Operation Manual Ver. 1.0.1

Firmware Ver. 6.0.2

Application Ver. 1.0.1.0

2020.12.28

Index

1	Overview	1
1.1	About this Operation Manual	1
1.2	Notes on using the application	1
1.2.1	About application	1
1.2.2	About connection	1
1.2.3	About Installation	1
2	Features	2
3	Setup	2
3.1	Setting APP permissions	2
4	Start the program	4
4.1	How to start	4
4.2	Startup screen (Signup screen)	4
4.3	Main screen	5
5	Connection	6
5.1	Registration of device in Smartphone	6
5.2	Registration of device in APP	8
6	Memory	11
6.1	Registration of Memory	11
7	Keypad/Voice	16
7.1	How to use Keypad input	16
7.2	How to use Voice input	18
8	Protocol	19
8.1	Edit	19
8.2	Protocol File	22
8.3	Modification of Protocol File	26
9	Logging function	28
9.1	Log File	28
10	Complete the operation	29
11	Troubleshooting	29

©2020 Icomes Lab Co., Ltd.

Icomes Lab, logo mark of Icomes Lab, pipetty and pipetty Pro are registered trademarks of Icomes Lab. The name "Android", the Android logo, "Google Play" and the Google Play logo are registered trademarks of Google LLC. In the United States and other countries. Bluetooth Low Energy is a registered trademark of Bluetooth SIG Inc. in the United States and is used by Icomes Lab Inc. under license.

In this Operation Manual, it may be abbreviated as follows.

Bluetooth Low Energy is Bluetooth LE or BLE

Other company names and product names are trademarks or registered trademarks of their respective companies.

The application specifications are subject to change without notice.

pipetty-Pro and pipetty Smart have built-in wireless equipment that has received construction design certification as a wireless station for low power data communication systems based on the Radio Law. EYSHCN: 001-A10745

1 Overview

1.1 About this Operation Manual

Please read this Operation Manual before using this product for the first time.
 We recommend that print this Operation Manual and store it near the product.
 This manual is for pipetty Smart or pipetty Pro Android™ smartphone apps only.
 Compatible with ver.1. *.* software (application) and ver.PP-6. *.* , ver.PS-6. *.* software (device).

1.2 Notes on using the application

1.2.1 About application

[Note]

- The log file is output as the operation result and does not guarantee the dispensed volume. Depending on the usage conditions, there may be a difference between the log file and the actual dispensed volume.
- This application display does not have a tracking function. Therefore, we cannot guarantee that the pipetting portion will be wrong due to the work.
- We are not responsible for any trouble or damage caused by using this application.
- The screen images in this manual are examples of use, so the contents may differ depending on the actual input contents.

1.2.2 About connection

[Note]

- In connection, when used around a wireless LAN, other wireless devices, around devices that emit radio waves such as microwave ovens, in places with many obstacles, or in other environments with poor radio wave conditions, the connection may be interrupted frequently, or communication may be lost. The speed may be extremely reduced, or an error may occur.
- If you repeatedly charge the battery without using it up, a phenomenon called the “memory effect” occurs in which the discharge voltage drops, and the connection may be interrupted or the communication speed may drop even after charging.
- For other precautions regarding wireless communication, please read the pipetty series Operation Manual.

1.2.3 About Installation

The following system is required for installation.

OS	Android™8.1 ~
----	---------------

2 Features


The "pipetty Smart App" which is interlocking with a pen-type electric pipette with a wireless communication function supports reduce researcher' s burden and make efficient operation by the reduce time required for switching between difference volume pipettes, mode switching, dispensing volume switching and experiment log saving.

3 Setup


The application can be downloaded from the "Google Play ™" link on our website.
[URL : <https://www.icomes.co.jp/product/pipetty-smart/>]

3.1 Setting APP permissions

* It is necessary to turn on the application permission from "Settings" on the smartphone start menu, "Apps and notifications" ⇒ "pipettySmartApp.Android" ⇒ "Permissions".



Turn on "Storage", "Microphone", and "Location" in APP permission screen. After completing the settings, tap ● "Close" on the smartphone.



4 Start the program

4.1 How to start

(1) On smartphone
Tap the "pipetty Smart" icon.



4.2 Startup screen (Signup screen)

When you run pipetty SmartApp, the startup screen (Registration screen) shown below is displayed.
Proceed the User registration by inputting information from this startup screen.

User registration is optional, but please be aware that if you do not register the necessary information, we may not be able to provide the service you desire.

The image shows two side-by-side screenshots of the Pipetty Smart App startup screen. The left screenshot shows the app title "Pipetty Smart App" and a blue "Sign up" button. The right screenshot shows the same app title and a registration form with the following fields: Name (icomes), E-mail (icomes@co.jp), Company or Lab (development), and Address (iwate.japan). Below the form are two blue buttons: "Revise" and "Registration". The "Registration" button is highlighted with an orange border.


Tap "Registration" to complete registration and move to the "Main" screen.

* Icomes Lab Co., Ltd. handles personal information obtained from customers in compliance with the Personal Information Protection Law and in accordance with the personal information protection policy.
Regarding protection of personal information, tap "Menu" ⇒ "Terms of Service" to open "About protection of personal information" on Icomes Lab Website, so be sure to confirm the contents.

4.3 Main screen


Description

- When you tap Menu, "Registration information" The "Terms of Service" screen will be displayed.



Be sure to read the Icomes Lab website.
 [Caution] Please use a terminal in the online environment.
 URL:
<https://www.icomes.co.jp/about/privacypolicy/#protection>

Please be sure to read "About the protection of personal information".




5 Connection

5.1 Registration of device in Smartphone


In order to use pipetty SmartApp, need to register pipetty device in smartphone.

- ① Open "Settings" in the start menu of the smartphone.
- ② Open "Settings" in the start menu of the smartphone.




- ③ On the connected terminal (Bluetooth), tap "Pair with new device".

* At this time, turn on the pipetty device (the display screen of the pipetty is flashing for mode selection).



- ④ "Nordic_UART" is displayed as "Available device" on the pairing screen with the new device.
- ⑤ "Nordic_UART" Tap to switch to the connected terminal screen, and "Nordic_UART" is displayed in "Currently connected terminals". (Pairing completed)
- ⑥ Tap "Nordic_UART" to switch to the detailed screen of the terminal and "Connect" will be highlighted.



- ⑦ Complete the pairing of pipetty and smartphone by "Back" or "Close" on the smartphone.


5.2 Registration of device in APP

Register the paired device in pipetty Smart App.

Register the device (20 μ L, 250 μ L, 1000 μ L) and the application side with same volume.

【Registration procedure (Registration example: pipetty 20 μ L)】


- ① Tap the setting of 20 μ L display on the Main screen.
- ② If it is first setting, "Incorrect pipette information" is displayed, so tap "OK".
- ③ Tap the information field of the device to register.



④ Confirm the MAC address to register from Select a device and tap it.

⑤ Confirm the connection preparation (make sure pipetty device is turned on), then tap "Yes".

⑥ Confirm that the MAC address is registered in the information field.





(Fig.1)

【How to confirm pipetty Number】

pipetty device Number can confirm by pressing the DOWN button while pressing the PUSH button while the mode selection is flashing.

The lower 6 digits of the MAC address are saved in the shipping state.

Also, can check the lower 6 digits of the factory-set MAC address by pressing the UP button while pressing the PUSH button.

Note: Cannot be displayed while connected to the application.


【About error messages】

● If the connection destination is incorrect

① When selecting a connection, if the device is incorrect, the following display will appear. Tap "OK" to reselect.

● If the pipetty is not turned on

(2) (3) The following display also appears when the pipetty is turned off. Turn on the pipetty and tap "OK" to reselect.



【Register multiple pipetty Smart device】


By registering up to 3 devices that are paired, can switch between capacities, modes, and dispensing volumes.

* Please note that cannot register multiple items with the same capacity device.

After registering 20 μ L, 250 μ L and 1000 μ L respectively, the screen will be as shown below.

① The button color will change when device registration is completed.

②③ When you tap the capacity button and the button turns blue, the operation is ready.




6 Memory



Up to 5 types of dispensing settings can be stored in memory.
By storing frequently used dispensing settings in memory, do not need to set every time.

6.1 How to register Memory

① Tap "Memory"



② To register from the Memory screen, tap "Edit" to enter the memory No. registration screen in ③.



- Memory No.1 ~ 5
- Select capacity
- Dispensing speed (5 steps)
- Aspiration speed (5 steps)
- Dispensing mode selects
Single-Mode: Dispense at once
Multi-Mode: Continuous dispensing
- Discharge volume setting
Single-Mode: Volume at once
Multi-Mode: Volume and number of times



* [Settings decimal point of dispensing volume]
 • 250 μ → 1.0-99.9 can be set. (However, the device display is up to 19.9)
 • 1000 μ → 5.0-99.9 can be set. (Same as above)
 Display on the device is up to 19.9, but the application can set and operate up to 99.9 after the decimal point.


④ Tap 2nd to register the dispensing setting in memory No.2

[Registration example: 20 μ L, Dispense speed 4, Aspiration speed 4, Multi-Mode, 1.5 μ L, 6 times]


⑤ Volume 20 μ L, Dispense speed 4, Aspiration speed 4, select ⑥ Multi-Mode




⑦⑧⑨ Tap the number of Multi-Mode dispensing volume μ L and the number of times Time to input (input example: Dispense 1.5 μ L 6 times)



⑩ After confirming the completion of input, tap "Return" ◀ to complete the dispensing setting registration and return to the Memory screen.



⑪ The registered dispensing volume is displayed, and tap the corresponding memory No.2, The setting volume will be transferred to the device.



(Fig.1)

* Other memory No. inputs can be set from "Edit" as well.




【Error message】

- Incorrect memory entry

①②③→When the input value of Volume is incorrect



④⑤⑥→When the input value is incorrect (Volume × number of times exceeds the specified value)



● In case of connection error

⑦ When the device is not turned on or when the main screen switch (blue) is not turned on



7 Keypad/Voice

Keypad input on the screen and voice input (Japanese / English) to set and enter the volume, dispensing mode, and dispensing volume of pipetty.

7.1 How to use Keypad

Pipetty Smart App

Menu

ON

20 μ L

250 μ L

1000 μ L

Memory

Keypad / Voice

Protocol

Pipetty Smart App

VOICE JP

20 μ L 250 μ L 1000 μ L

0.1 - 20 1 - 250 5 - 1000

Single-mode Multi-mode

0 μ L 0 μ L 0^{time}

1 2 3

4 5 6

7 8 9

. 0 BS

Send

Voice input ON and JP / EN switching

Select capacity of connected devices

Select dispensing mode
Single-Mode: Dispense at once
Multi-Mode: Continuous dispensing

Dispensing volume setting
Single-Mode: Dispense volume at 1 time
Multi-Mode: Dispense volume and number of times

Keypad (numeric input)

Enter the setting

[Keypad input registration example: 1000 μ L, Single-Mode, 500 μ L]

- ① Tap 1000 μ L capacity and select Single-Mode in dispensing mode
- ② When tap the number part of the dispensing volume setting, the number disappear, and can proceed input.

① Pipetty Smart App

VOICE JP

20 μ L 250 μ L 1000 μ L

0.1 - 20 1 - 250 5 - 1000

Single-mode Multi-mode

0 μ L 0 μ L 0^{time}

1 2 3

4 5 6

7 8 9

. 0 BS

Send

② Pipetty Smart App

VOICE JP

20 μ L 250 μ L 1000 μ L

0.1 - 20 1 - 250 5 - 1000

Single-mode Multi-mode

0 μ L 0 μ L 0^{time}

1 2 3

4 5 6

7 8 9

. 0 BS


Send

* [Settings below the decimal point of dispensing volume]

- 250 μ → 1.0 - 99.9 can be set. (However, the display of device is up to 19.9)
- 1000 μ → 5.0 - 99.9 can be set. (Same as above)

The display on device is up to 19.9, but the application can set and operate up to 99.9 after the decimal point.

③ After inputting a numerical value on the keypad, tap "Send" to transfer to device.




* Tap the "Back button" on the smartphone to reset the entered information.

[Error messages]


- If the keypad value is entered incorrectly

① When the input value of Volume is incorrect

② When Volume × number of times exceeds the specified value



Input the "Dispense volume" which does not exceed the specified value.



Input the "Dispense volume × Number of times" which does not exceed the specified value.

7.2 How to use Voice input

[Voice input (Japanese)]


① Turn on the capacity button from the Main screen

② Turn on the VOICE button and tap "JP (Japanese)".

[Input example: Capacity 1000 μ L, Multi-mode, 100 μ L, when inputting voice for setting 5 times (Japanese)]

③ Speak "1000 micro" into the microphone of your smartphone.

→ The recognized words are displayed as a keypad and the capacity switches to 1000 μ L.



④ When you speak "multi", the keypad is displayed, and it switches to Multi-mode.

⑤ The setting of the dispensing volume is called "volume".

⑥ When the value is blank, speak the value "100". When recognized, the dispensing volume in μ L is displayed.



- ⑦ It can also specify the number of dispenses in the same way.
- ⑧ The numerical value is called "5".
- ⑨ When you say "Send" after completed the input, the data will be transferred to device.




8 Protocol


8.1 Edit

Create and input protocol files and edit.


- ① Tap "Edit" to open the edit screen.
- ② If there is no target file for "Save File" and "Delete File", an error will be displayed.



③ - ⑤ Tap "Add" to add items and create the required protocol.



⑥ - ⑧ Select "Capacity" and "Mode"



* [Select Mode]

" Single" → Quantitative dispensing operation "Multi" → Equal volume continuous dispensing operation

"Super Multi" → Different volume continuous dispensing operation


" Mix" → Mixing (Pipetting) operation "None" → No selection (When do not run device)

* [Selecting the number of times]

Only "Multi" can select the number of times.

"Mix" repeats aspiration and dispensing while pressing the operation button of device. (One-time operation)

⑨ - ⑪ Then enter the "Volume (Dispensing volume)". Single mode is automatically set to once.



* [Settings below the decimal point of dispensing volume]


- 250 μ → 1.0 - 99.9 can be set. (However, the display of device is up to 19.9)
- 1000 μ → 5.0 - 99.9 can be set. (Same as above)

The display on device is up to 19.9, but the application can set and operate up to 99.9 after the decimal point.


⑫ Continue tapping "Add" to add necessary items and input.

(* In the case of Super Multi mode, enter the volume to be dispensed into each item once.)

⑬ After entered the file name, "Save File" is enabled.




- ⑮ Tap "Save File" to complete and save to a file.
- ⑯ If the same file name exists, a confirmation screen for overwriting will be displayed.
 OK ⇒ Overwritten
 Cancel ⇒ Please enter a new name




8.2 File (Call protocol file)

Calls the entered protocol file and transfers it to device. (Example: Protocol 1)


- ⑰ Before opening File, turn on the connection button.




② Tap File first. And tap ③ "File Name" twice to display ④ "Please select a file".



⑤ Select the specified file and tap "Start" to start the protocol.
Press the PUSH button on device, and after operation in item 1 is complete, check mark will appear in "Done".




- ⑧ If enter only the protocol characters as in No. 4, tap "Done" manually.
 - ⑨ A check mark is added and be able to proceed to the next section.
 - ⑩ When the last check mark is completed, a PDF file will be output as a record in the log file in the smartphone.
- * If operation canceled in the middle of the protocol, a PDF file that records the time until the check mark is output.




[Display of popup]

- ① If the following message is displayed during protocol operation, press the "PUSH" button on device to turn on.
- ② When Cancel is selected, ③ "Failed to connect device" is displayed and the screen returns to the protocol screen. If so, start the protocol over again.




④ If the following message appears while the protocol is operating, and the message ⑤ "... mode start failed" appears after tapping "Yes", the device may not be in the mode selection state. Set the device to standby state and try the protocol again.



[In case communication is interrupted]


- If "Incorrect pipette information" or "Failed to connect" is displayed on the application screen, please try again from Section 5.2 (pipetty Smart registration). In that case, please remove the battery from device once, and re-insert it.



8.3 Modify the protocol file




When need to modify the item of the called file, tap Edit and "Tap File Name" twice.

*If tap only once, the screen "Enter file name" will appear, and the protocol file will be created. Return with Cancel and tap twice to try again.




[File modification]


- ① Tap Edit.
- ② Tap the File Name field twice.
- ③ Tap the file name need to modify.



- ④ The selected protocol is displayed.
- ⑤ Tap twice on the part need to modify. (Example: Volume in Section 7 changed from 160 to 150 μL)
- ⑥ Enter the value to be modified.



- ⑦ Confirm that the value has been changed, and tap "Save File".
- ⑧ A confirmation message for overwriting the file will be displayed. OK \Rightarrow Overwrite.
If set a new file name, Cancel \Rightarrow Return, tap the File Name field once to enter a new name, and tap OK (⑨ is the correction completion screen)



9 Logging function

In the logging function mode, basic operations are performed by device, and only the operation results are recorded in the log.

9.1 Log files

Operation result log is from the shared storage in the smartphone file.

The execution result log is output as a pdf file to the shared storage

\ Android \ data \ jp.co.icomes.pipettysmartapp \ files \ PipettySmartApp

in the smartphone file. The following is the output of the operation result by the protocol of 8.2.

Protocol A

No	Protocol	Capacity	Mode	Volume	Times	EndTime
1	Dispense 500 micro of Reagent A into the tube	1000µ L	Single	500.0	1	2020/06/03 11:17:37
2	Dispense 20 micro of Reagent B	20µ L	Single	20.0	1	2020/06/03 11:17:51
3	Dispense 5 micro of mixed liquid into wells 4times	20µ L	Multi	5.0	4	2020/06/03 11:18:08
4	Leave for 1 minute	-	-	-	-	2020/06/03 11:18:09
5	Dispense Reagent C into wells by 100,125,150, 175 micro	1000µ L	SuperMulti	100.0	1	2020/06/03 11:18:15
6		1000µ L	SuperMulti	125.0	1	2020/06/03 11:18:17
7		1000µ L	SuperMulti	160.0	1	2020/06/03 11:18:19
8		1000µ L	SuperMulti	175.0	1	2020/06/03 11:18:26

The following operations are saved as the operation result log.

Protocol (Item) Capacity (Device volume) Mode Volume (Dispensing volume)

Times (Number of times) End Time (Completion recording time)

* EndTime is output for the part with the "Done" check mark on protocol operation.

[File storage destination for protocol and memory registration]

Shared storage \ Android \ data \ jp.co.icomes.pipettysmartapp \ files \ PipettySmartApp.

【Notes】


Depending on the smartphone model, the "Protocol file", "Memory file", and "Log file" may not be viewable from the PC. In that case, please restart your smartphone or tablet device.

【Notes】

If an error (Err004) is displayed on device during aspiration, the device has exceeded the normal operating temperature range, so suspend the operation and place the device in a cool place for a certain period of time or change the dispensing volume.

10 Complete the operation

After completed the operation, return to main screen and tap button of used capacity to disconnection.



11 Troubleshooting

If the problem cannot be solved by the following method or if any other abnormality occurs, please contact us or distributor.

Problem	Inferred factors	Solutions
Cannot connect to device	No battery inserted in device	Insert the battery
	The device display is not the mode setting screen	Switch the display of device to the mode setting screen (initial screen)
	Low battery	Replace with a charged battery
	Battery has memory effect	After discharging the battery on the 1st scale as a guide, repeat charging 2-3 times.
	Battery has reached cycle life	Replace with a new battery
	Battery is inserted in the opposite direction	Insert the battery correctly
	No location permissions in Smartphone	Setting location information authority of smartphone
Easy to disconnect	Low battery	Replace with a charged battery
	Battery has memory effect	After discharging the battery on the 1st scale as a guide, repeat charging 2-3 times.
	Battery has reach	Replace with a new battery
	Affected by peripherals	Check if there is any influence of other wireless communication

Problem	Inferred factors	Solutions
Easy to disconnect	The device and smartphone are too far apart	Reduce the distance between device and the smartphone
	There is a wall or other shield between device and the smartphone	Make sure there is no obstruction between device and smartphone
	The timing to operate is early	Operate according to the display of the application
Cannot find protocol data / log data	Storage settings are hidden in Smartphone	Show storage settings of Smartphone
Cannot input voice	No microphone authority	Set microphone permission on Smartphone
Cannot create folder	There is no permission for storage (Create log)	Set storage permissions on Smartphone
"Incorrect pipette information" is displayed	Incorrect settings or device selection	Check the device capacity and MAC address and try again in Section 5.2.
"Connection failed" is displayed	Device is not turned on	Press the PUSH button to turn on
	Low battery	Replace with a charged battery
	Problems on the smartphone side	Retry the "pipetty smartphone registration" in section 5.1.

Customer Support



Icomes Lab Co., Ltd.

2-4-23, Kitaiioka, Morioka, Iwate, 020-0857 JAPAN

Healthtech Innovation Hub

TEL: +81-19-601-8228 FAX: +81-19-601-8227

E-mail: globalinfo@icomes.co.jp Website: <https://www.icomes.co.jp/en>