

Handy type **Tofutty** Electric Dispenser

**PATENT
PENDING**

Features

Compact and Cordless body

- Easy movement in narrow space
- Small size (height 149mm , Width 22.5mm , depth 38mm)
- Light Weight (75g)

Precise and Easy operation

- $\Phi 8$ micro actuator technology of icomes makes it possible to apply low liquid of high viscosity materials in high precision.

Wireless Application (Optional function)

- Wireless application automatically makes a record of your dispensing work.

Shake-Prevention During dispensing (Optional function)

- It can be connected between dispenser and foot switch by wireless communication.





Product Price ¥75,000

Syringe Price (20set) ¥12,000

Specification

	Specification
Size	H: 149mm , W: 22.5mm , D: 38mm
Weight	75g
Capacity	1,000 μ L(1cc)
Power	AAA battery (Ni-MH) (Specify : Panasonic rechargeable battery EVOLTA)
Syringe nozzle	General-purpose needles are available (Recommendation : Musashi Engineering taper nozzle 22G)
Viscosity	Under 20,000mPa
Repeatability(CV)	Coefficient of variation:10%(1 μ L)
Syringe	Clear syringe (Material:Polypropylene) UV block syringe (Material:Polypropylene)

Accuracy of liquid discharge




Option

Wireless communication



Wireless communication automatically record date time and dispense volume.

Foot switch



Foot switch prevents hand's shaking control at the time the button is pushed.



Location : Kitaiioka 1-8-25, Morioka-shi,
Iwate 020-0857, Japan

TEL : +81-19-601-8228 FAX : +81-19-601-8227

E-mail : info@icomes.co.jp

URL : <http://www.icomes.co.jp>